

## Actril® Cold Sterilants

### Research Report: Clostridium Difficile Endospores and PAA Germicides

#### Introduction

Actril Cold Sterilant is a Peracetic Acid (PAA) based sterilant that has specific claims against even the hardest class of organisms – spores. While Actril Cold Sterilant has been shown effective against the AOAC spore test organism – *Bacillus subtilis* – this white paper examines the effectiveness of current germicides, as well as, Actril Cold Sterilant against *Clostridium difficile*.

#### History Of *C. difficile*

While identified in the literature for over 70 years, *C. difficile* began to be increasingly recognized in the 1970s as an organism that had adapted and was becoming progressively more antibiotic resistant. *C. difficile* is a spore forming, gram positive bacteria that is highly resistant to acidic environments such as seen in the gastrointestinal (GI) system. Because they can thrive in the intestinal tract, *C. difficile* releases toxins directly into the GI system, destroying the intestinal lining resulting in diarrhea and extreme dehydration.<sup>1</sup>

From 1999 to 2004, *C. difficile* related deaths almost quintupled with the death rate increasing by 35% per year in the United States.<sup>2</sup> In other nations, *C. difficile* has surpassed MRSA as the “superbug” of greatest concern in terms of numbers of deaths and continued growth rates.<sup>3</sup>

#### Germicides and *C. difficile*

*C. difficile* has two states: vegetative and endospore. In the vegetative state, *C. difficile* is readily destroyed by a large number of germicides such as quaternary ammoniums, dilute hypochlorite (bleach) solutions (1%), phenols, and alcohols.<sup>4</sup> In contrast, *C. difficile* in the endospore state is characterized by a thickened cell wall, which provides protection against drying-out and acidic environments. This cellular structure enables the *C. difficile* endospore to be resistant to the previously-mentioned germicides. In a study reported by Dr. William Rutala in 2006, the following disinfectants demonstrated no measurable activity at

20 minutes against *C. difficile* spores:

- (1) Chlorhexadine,
- (2) Vesphene (phenol),
- (3) 70% isopropyl alcohol,
- (4) 95% ethanol,
- (5) 3% hydrogen peroxide,
- (6) Clorox disinfecting spray (65% ethanol, 0.6% quaternary ammonium),
- (7) Novaplus (10% povidone iodine) and
- (8) Virox’s Accel (0.5% hydrogen peroxide).<sup>5</sup>

In fact, in one study performed in England, two disinfectant compounds – quaternary ammonium and hydrogen peroxide – appeared to even encourage the growth of *C. difficile* spores.<sup>6</sup>

#### Actril vs. Spores

Actril Cold Sterilant is a broad spectrum germicide with sporicidal claims. In testing conducted under AOAC Sporicidal Test Protocols, Actril Cold Sterilant was successful with complete kill against *Bacillus subtilis* and *Clostridium sporogenes*.<sup>8</sup>

Actril Cold Sterilant has been shown to be effective against spores in general, a more specific test was conducted to look at its effectiveness in a short time period against *C. difficile* spores.<sup>9</sup> This test is summarized below:

Test Methodology: ASTM E 2197-02 – Conducted with ATCC 700792 (*C. difficile* spore)

Test Results: The inoculated carriers had a starting population of  $4 \times 10^5$  cfu. 10 carriers were exposed for 10 minutes to Actril solution. The results at the end of the ten minute period were:

Germicide	Carriers	% Reduction
Actril Cold Sterilant	10	99.999

## Discussion

While a few antibiotic regimes such as vancomycin continue to be effective for the most part against *C. difficile* spores, good infection control practices dictate environmental reduction of the organism before infection. As previously indicated, not all germicides are effective against *C. difficile* spores and some may even promote growth.

In testing by Minntech Corporation, Actril Cold Sterilant significantly reduced the populations of *C. difficile* spores 99.999% and 99.9%, respectively.

## References

- 1 McMaster-Baxter Nicole Pharm.D., Musher Daniel M. MD. Clostridium difficile: Recent Epidemiologic Findings and Advances in Therapy. Pharmacotherapy. 2007;27(7):1029-1039.
- 2 Redeling Matthew D., Sorvillo Frank, Mascola Laurene. Increase in Clostridium difficile-related Mortality Rates, United States, 1999-2004. Emerging Infectious Disease. 2007; 13(9):1417-1419.
- 3 [www.statistics.gov.uk/pdffdir/deaths0207.pdf](http://www.statistics.gov.uk/pdffdir/deaths0207.pdf)
- 4 From EPA labels.
- 5 Rutula William, et. al. Surface Disinfection: New Processes and Products. Presentation at 2006 APIC National Meeting. 2006.
- 6 Fawley Warren N., et.al. Efficacy of Hospital Cleaning Agents and Germicides Against Epidemic Clostridium difficile Strains. Infect Control Hosp Epidemiol. 2007 Aug ;28 (8):920-5.
- 7 Minncare Cold Sterilant: Research Data Report. 2003.
- 8 Actril Cold Sterilant: Research Data Report. 2003.
- 9 Maltais Jo-Ann, Stern Thomas. Technical Report: A Superior Biocide for Disinfecting Reverse Osmosis Systems. 2003.
- 10 Minntech Study. C.difficile and Actril Cold Sterilant. May 2008.

Actril® is a registered trademark of Minntech Corporation, a Cantel Medical Company



Mar Cor Purification 4450 Township Line Road Skipack, PA 19474-1429 Tel: (484) 991-0220 Toll Free: (800) 346-0365 Fax: (484) 991-0230	Mar Cor Purification 14550 28th Avenue North Plymouth, MN 55447 Tel: (484) 991-0220 Toll Free: (800) 633-3080 Fax: (763) 210-3868	Mar Cor Purification 3250 Harvester Road - Unit 6 Burlington, ON L7N 3W9 Tel: (905) 639-7025 Toll Free: (800) 268-5035 Fax: (905) 639-0425	Mar Cor Purification Sourethweg 11 6422 PC Heerlen The Netherlands Tel: (+31) 45 5471 471 Fax: (+31) 45 5429 695	Mar Cor Purification 1A International Business Park, #05-01 Singapore 609933 Tel: (+65) 6227 9698 Fax: (+65) 6225 6848
--	--	---	---	---